

SPENGLER MASTER PALM 2 PORTABLE OXIMETER + SILICONE PROTECTIONS

Ref. 223000P



18/08/2026 14:01:22 – For information purposes only, can be subject to change

SPENGLER MASTER PALM 2 PORTABLE OXIMETER + SILICONE PROTECTIONS

Ref. 223000P



Spengler Master Palm 2 pulse oximeter: reliable SpO₂ monitoring with silicone covers for enhanced durability.

Description

Spengler Master Palm 2 pulse oximeter: reliable SpO₂ monitoring with silicone covers for enhanced durability.

Features

CE Marking	CE	Medical device	Class IIb
Bar code	3661040222608	HS Code	9018199000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.718
Bare Unit / Weight (kg)	0.105		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES