

PAIR OF SCREW TYPE EARTIPS/BLACK/FOR "TRIANGLE" STETHOSCOPE

Ref. 506920



Ear tips, designed for safe professional use.

Description

Spare part (x2) for Cardio Prestige® Cardio Sonic® Pédiasonic® Magister® Laubry® Pulse® & Dual Pulse® stethoscopes

18/08/2026 13:26:36 – For information purposes only, can be subject to change

PAIR OF SCREW TYPE EARTIPS/BLACK/FOR "TRIANGLE" STETHOSCOPE



Ref. 506920

Features

CE Marking	NONE	Medical device	No
Bar code	3700446009410	HS Code	90330000
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		