



Ref. 507219



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
5 ANS
5 YEARS
GUARANTEE



18/08/2026 13:32:34 – For information purposes only, can be subject to change

=====

Ref. 507219

Spengler takes care of newborns!

Description

The dual Magister Neonatal est dedicated to newborn auscultation.

With its dual head chestpiece equipped of a floating diaphragm on the large side, Magister Neonatal is perfectly adapted to babies for listening high and low frequencies.

Its PVC non-chill rings on both sides of the chestpiece allows a perfect comfort for the newborns.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446042325	HS Code	90181990
Warranty (months)	60		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Bare Unit / Length (cm)	76
Packed unit / Weight	0.25		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	NA	Maximum humidity	95

OTHER FEATURES