

## DUAL PULSE STETHOSCOPE / BLACK

Ref. 507400



## DUAL PULSE STETHOSCOPE / BLACK

Ref. 507400



### Description

Pearly tube model. For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaural with built-in spring.

### Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | CE            | Medical device            | Class I  |
| Bar code                           | 3700446009199 | HS Code                   | 90181990 |
| Warranty (months)                  | 36            |                           |          |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | 11            | Packed unit / Height (cm) | 4        |
| Packed unit / Length (cm)          | 24            | Packed unit / Weight      | 0.237    |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | NA            | Max temperature           | NA       |
| Minimum humidity                   | NA            | Maximum humidity          | NA       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Colors                             | Black         | Packed unit / Quantity    | 1        |