

DIAPHRAGM/BLACK RETAINING RING FOR LAUBRY

Ref. 507936



Description

Spare part for Laubry® stethoscopes

18/08/2026 12:20:02 – For information purposes only, can be subject to change

**DIAPHRAGM/BLACK RETAINING RING FOR
LAUBRY****Ref. 507936****Features**

CE Marking	NONE	Medical device	Class I
Bar code	3700446009670	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	4	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	4	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		