

VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)

Ref. 512220



CE
0459



Description

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

18/08/2026 23:37:23 – For information purposes only, can be subject to change

**VAQUEZ LAUBRY NANO / VELCRO CUFF /
BLACK / CHILD (S)****Ref. 512220****Features**

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007256	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.582
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		