

4VS MONITOR NIBP

Ref. 551600



CE
0123



18/08/2026 12:19:27 – For information purposes only, can be subject to change

4VS MONITOR NIBP

Ref. 551600



Practical for hospitals thanks to its base it is easy to move it or even position it near the patient's bed

Description

The 4VS PNI vital station has an LCD screen for better reading, visual and audible alarms are also triggered when the parameter has exceeded its limit.

It offers you remarkable reliability and accuracy for nearly 10 hours, it is possible to recharge it on mains. •

1 Vital signs: PNI•

3 Blood pressure measurement modes, infants, child and adult•

Possibility to program physiological

alarms• 4 sizes of armbands: neonatal, pediatric, adult, large adult (delivered with adult cuff)

• 3500 groups of data saved•

Battery: 10h• Size: 180 mm .

250 mm . 180 mm; Weight: 2.0 kg (with battery)

• Roller stand, wall mount, optional printer

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446032944	HS Code	9018199560
Warranty (months)	18	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.76
Bare Unit / Weight (kg)	2		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--