

NIPB TUBING 2,5M FOR 4VS MONITOR

Ref. 551624



Features

CE Marking	NONE	Medical device	No
Bar code	3661040213842	HS Code	85049099
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	16	Packed unit / Height (cm)	3
Packed unit / Length (cm)	17	Packed unit / Weight	0.118

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---