


SPO2 ADULT SENSOR FOR 4VS MONITOR

Ref. 551630



CE
0537



18/08/2026 11:55:26 – For information purposes only, can be subject to change

SPO2 ADULT SENSOR FOR 4VS MONITOR

Ref. 551630



Rigid adult SpO₂ sensor for Spengler 4VS, designed for reliable and continuous monitoring.

Description

SPO2Adult Rigid Sensor for SPENGLER 4 VS Vital Signs Monitor

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031176	HS Code	90181990
Warranty (months)	6	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.103

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---