


SPO2 SENSOR WRAP FOR 4VS MONITOR

Ref. 551633



CE
0537



Flexible and reliable sensor for comfortable SpO₂ monitoring with the 4VS monitor.

Description

SPO2 Sensor Wrap for SPENGLER 4 VS Vital Signs Monitor

18/08/2026 11:53:56 – For information purposes only, can be subject to change

SPO2 SENSOR WRAP FOR 4VS MONITOR

Ref. 551633



Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031206	HS Code	90181990
Country	Uruguay		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.084

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---