

SPO2 PEDIATRIC SILICON SENSOR FOR 4VS MONITOR

Ref. 551635



CE
0537



A soft and accurate sensor for pediatric oxygen monitoring with the Spengler 4VS monitor.

Description

SPENGLER Pediatric Silicone SPO2 Sensor for Vital Signs Monitor 4 VS

18/08/2026 13:55:14 – For information purposes only, can be subject to change

SPO2 PEDIATRIC SILICON SENSOR FOR 4VS MONITOR

Ref. 551635



Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031220	HS Code	90181990
Country	Uruguay		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.083
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		