


SPO2 EAR CLIP SENSOR FOR 4VS MONITOR

Ref. 551636



CE
0537



Description

SPO2 ear sensor for SPENGLER 4 VS vital signs monitor

18/08/2026 12:41:22 – For information purposes only, can be subject to change

SPO2 EAR CLIP SENSOR FOR 4VS MONITOR

Ref. 551636



Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031244	HS Code	90181990
Country	Uruguay		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.082
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		