

CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH

Ref. 570150



18/08/2026 12:40:39 – For information purposes only, can be subject to change

CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH

Ref. 570150



Enhanced ENT examination with the Spengler Clar mirror, supplied with a protective carrying case.

Description

Suitable for all head sizes thanks to a flexible metal headband.
 Indirect lighting by reflecting mirror.
 Suited to lighting required in daily practice.
 Built-in battery compartment, two CR 123 batteries
 Supplied with carrying bag

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446011390	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Packed unit / Height (cm)	11
Packed unit / Length (cm)	24	Packed unit / Weight	0.655
Bare Unit / Weight (kg)	0.23		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		