

CABLE USB-C - USB MALE 2M POUR SMARTLED 5500-R



Ref. 579535

A sturdy and convenient cable for easy connection to the SMARTLED 5500-R.

Description

This 2-meter USB-C to USB male cable ensures fast and reliable connectivity with the SMARTLED 5500-R monitor. Its extended length provides flexibility in setup while maintaining stable data transmission.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040216782	HS Code	9018110000
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	17	Packed unit / Weight	354
Bare Unit / Weight (kg)	354		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	80
Minimum humidity	5	Maximum humidity	90

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---