

PROBE FOR FETAL DOPPLER / 2MHZ

Ref. 743310



CE
0123



Spengler 2 MHz probe for fetal Doppler, ideal for listening to the baby's heartbeat.

Description

2 MHz probe - fetal probe for doppler

18/08/2026 12:36:22 – For information purposes only, can be subject to change

PROBE FOR FETAL DOPPLER / 2MHZ

Ref. 743310



Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446008185	HS Code	90181990
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	3	Packed unit / Height (cm)	2
Packed unit / Length (cm)	13	Packed unit / Weight	0.145

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	55
Minimum humidity	5	Maximum humidity	93

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---