

PROBE FOR VASCULAR DOPPLER / 5 MHZ

Ref. 743311



CE
0123



A precise and ergonomic 5 MHz probe for vascular assessments with the Spengler Doppler.

Description

5 MHz probe - vascular probe for doppler

18/08/2026 13:36:54 – For information purposes only, can be subject to change

PROBE FOR VASCULAR DOPPLER / 5 MHZ**Ref. 743311****Features**

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446008178	HS Code	90181990
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	3	Packed unit / Height (cm)	2
Packed unit / Length (cm)	13	Packed unit / Weight	0.145

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	55
Minimum humidity	5	Maximum humidity	93

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---