

PROBE FOR VASCULAR DOPPLER / 8 MHZ

Ref. 743312



CE
0123



High-precision 8 MHz probe, suitable for cardiac and vascular examinations with Spengler Doppler.

Description

8 MHz probe - cardiac probe for doppler

18/08/2026 13:54:54 – For information purposes only, can be subject to change

PROBE FOR VASCULAR DOPPLER / 8 MHZ**Ref. 743312****Features**

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446008161	HS Code	90181990
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	3	Packed unit / Height (cm)	2
Packed unit / Length (cm)	13	Packed unit / Weight	0.145

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	55
Minimum humidity	5	Maximum humidity	93

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---