

LANCET 0,32MM 30G (200 UNITS)

Ref. 745013



CE
0123



LANCET 0,32MM 30G (200 UNITS)**Ref. 745013**

=====

Description

Lances for use with the TB 100 Lancing Pen.

-

Needle diameter 30G = 0.298mm ~ 0.320mm

-Sterile, single use - Length 32mm +/- 1mm- Width 6.5mm +/-

0.5mm

-Outer needle 3.3mm +/- 0.3mm- Packaging: box of 200 (

20 sachets of 10 lancets)

Find our complete kit Glucometer SPENGLER TB 100 (ref. 745000)

Refund LPP code: 6111373

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041687	HS Code	90183190
Warranty (months)	24	Sterile	Yes

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	5
Packed unit / Length (cm)	10	Packed unit / Weight	0.22

SPECIFIC STORAGE CONDITIONS

OTHER FEATURES

Code LPP	6111373	Packed unit / Quantity	1
----------	---------	------------------------	---