


LUXASCOPE AURIS LED 2.5 V - LUXAMED

Ref. A1.416.214

LUXAMED®



CE



LUXASCOPE AURIS LED 2.5 V - LUXAMED

Ref. A1.416.214

LUXAMED®

Innovative otoscope in pocket size

Description

Illumination: LED-Ring with 6 LEDs (concentric),
no fibre optic, field of vision is not restricted

- Maintenance-free, because no lamp replacement necessary
- Can be dimmed
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours
by constant current drive
- Intensity of illumination: approx. 10.000 Lux
- Colour temperature: approx. 4.200 K
- Optimized design of head
- Pivoting 3-fold magnifying glass with locking function
- Insufflation port for pneumatic test
- Two-component handle (aluminium / fibre-glass reinforced),
also combinable with LuxaScope LED Dermatoscope- and
Ophthalmoscope head
- Electronic ON/OFF push-button switch
- Power supply: 2x AAA Alkaline batteries, inclusive
- Incl. 20 disposable ear funnels, grey,
each 10 x 2.5 mm and 10 x 4.0 mm

LUXASCOPE AURIS LED 2.5 V - LUXAMEDRef. A1.416.214**LUXAMED®****Features**

CE Marking	CE	Medical device	Class I
Bar code	04260381960027	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	11	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.12		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Blue
--------	------