

OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

Ref. A1.426.114

LUXAMED®



Rechargeable USB



OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

LUXAMED®Ref. A1.426.114

Innovative otoscope in pocket size

Description

Illumination: LED-Ring with 6 LEDs (concentric),
no fibre optic, field of vision is not restricted

- Maintenance-free, because no lamp replacement necessary
- Can be dimmed
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours
by constant current drive
- Intensity of illumination: approx. 10.000 Lux
- Colour temperature: approx. 4.200 K
- Optimized design of head
- Pivoting 3-fold magnifying glass with locking function
- Insufflation port for pneumatic test
- Two-component handle (aluminium / fibre-glass reinforced),
also combinable with LuxaScope LED dermatoscope- and
ophthalmoscope head
- Electronic ON/OFF push-button switch
- Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive
- Incl. 20 disposable ear funnels, grey, each 10 x 2.5 mm and 10 x 4.0 mm
- Incl. recharger with USB plug (adapter EU/UK/US)

Available in :

- White : A1.426.914
- Blue : A1.426.214
- Grey : A1.426.314
- Black : A1.426.114

OTOSCOPE LUXASCOPE AURIS LED 3.7 V - LUXAMED

LUXAMED®
Ref. A1.426.114

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960041	HS Code	90189084
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.45
Bare Unit / Weight (kg)	0.105		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		