


LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

Ref. A1.616.114

LUXAMED®



Technologie CCT
2 couleurs de lumière



18/08/2026 11:51:24 – For information purposes only, can be subject to change

LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

LUXAMED®Ref. A1.616.114

Innovative otoscope in pocket size

Description

Innovative illumination: LED-Ring with 2x4 LEDs (concentric)
switchable with one button press : warm white light or cold white light,
no fibre optic,
field of vision is not restricted

- Maintenance-free, because no lamp replacement necessary
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours due to ideal operation by constant current drive
- Colour temperature: ring 1 approx. 3.000 K, ring 2 approx. 6.000 K
- Optimized design of head
- Pivoting 3-fold magnifying glass with locking function
- Insufflation port for pneumatic test
- Two-component handle (aluminium/fibre-glass reinforced), also combinable with Luxascope Dermatoscope CCT LED head
- Electronic ON/OFF push-button switch
- Power supply: 2x AAA Alkaline batteries, inclusive
- Incl. 20 disposable ear funnels, grey, each 10 x 2.5 mm and 10 x 4.0 mm

Available in :

Black : A1.616.114

Blue : A1.616.214

Grey : A1.616.314

White : A1.616.914

LUXASCOPE AURIS CCT LED 2.5 V - LUXAMED

LUXAMED®
Ref. A1.616.114

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963578	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		