

**LUXASCOPE AURIS SET LED 2.5 V  
"COLOUR-EDITION", MAGENTA  
(OTOSCOPE + OPHTHALMOSCOPE)**

Ref. AF.416.974

**LUXAMED®**



CE



Complete set for precision diagnosis: Auris LED 2.5V Magenta otoscope + LuxaScope LED 2.5V Magenta ophthalmoscope

## Description

18/08/2026 13:57:35 – For information purposes only, can be subject to change

## **LUXASCOPE AURIS SET LED 2.5 V "COLOUR-EDITION", MAGENTA (OTOSCOPE + OPHTHALMOSCOPE)**

**LUXAMED®**

Ref. AF.416.974

Otoscope :

Lighting: LED ring with 6 LEDs (concentric),

no fibre optics, the field of vision is not limited.

- Maintenance-free, as no bulb change is required.
- Possibility of varying brightness intensity.
- Automatic switch-off after 3 minutes
- LED life: approx. 100,000 hours at constant current

current

- Illuminance: approx. 10,000 Lux

- Colour temperature: approx. 4,200 K

- Optimised head design

- Swivelling magnifier, 3x magnification

- Injection port for pneumatic test

- Handle made from two materials: aluminium and fibreglass

fibreglass, compatible with the head of the Lux dermatoscope and

LuxaScope ophthalmoscope head

- ON/OFF push-button switch

- Power supply: 2x AAA alkaline batteries, included

- Includes 20 grey disposable specula, 10 of which are 2.5mm and 10 of which are 4.0mm specula

Ophthalmoscope :

- Correction values: + 20 bis - 20 D ( $\pm$  1 2 3 4 6 8 10 15 20 D)

- 5 different filters available: large spot, small spot, green filter, half-moon and reticle

- Illumination: LED

# LUXASCOPE AURIS SET LED 2.5 V "COLOUR-EDITION", MAGENTA (OTOSCOPE + OPHTHALMOSCOPE)

**LUXAMED®**

**Ref. AF.416.974**

- Brightness can be dimmed to avoid glare glare
  - Automatic switch-off after 3 minutes
  - LED life: approx. 50,000 hours at constant current current
  - Colour temperature: approx. 4,200 K
  - Optimised head design
  - Handle made from two materials: aluminium and fibreglass reinforced, compatible with otoscope and dermatoscope heads
- LuxaScope LED
- ON/OFF push-button switch
  - Includes 2 AAA alkaline batteries

## Features

|            |                |                   |         |
|------------|----------------|-------------------|---------|
| CE Marking | CE             | Medical device    | Class I |
| Bar code   | 04260381964490 | Warranty (months) | 24      |
| Country    | Germany        |                   |         |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -10 | Max temperature  | 45 |
| Minimum humidity    | 10  | Maximum humidity | 90 |

### OTHER FEATURES