


EMERGENCY/TRANSPORT MONITOR C30 STANDARD CONFIGURATION COMEN

Ref. BCMONI07

COMEN



CE
1639



DISPOSITIF MEDICAL
CLASSE IIb
MEDICAL DEVICE

GARANTIE
3
3 YEARS
GUARANTEE



18/08/2026 13:37:33 – For information purposes only, can be subject to change

EMERGENCY/TRANSPORT MONITOR C30 STANDARD CONFIGURATION COMEN

COMEN

Ref. BCMONI07

A 5-in-1 monitor with monitoring of 5 vital parameters. Intuitive thanks to its color touch screen and very easy to carry because ultra compact.

Description

The C30 multi-parameter monitor, dedicated to first aid for on-site monitoring or in emergency vehicles. It measures non-invasive blood pressure (PNI), oxygen saturation (SpO2), respiratory frequency (RESP), capnography (EtCO2), body temperature and has a 3-lead ECG.

Parameters: PNI, SpO2, ECG, Respiration, 2 temperature channels (probes to be added), EtCO2 (sensor to be added).

Screen: 4.3" touch with backlight.

Double alarm: audible (adjustable) and 3-level visual with 3-color smart alarm system (blue, yellow, red).

PNI: suitable for all types of patients (adult, child and newborn).

Storage: 48 hours of data storage and trend analysis.

Battery: Built-in 2200mAh lithium (4 hours battery life).

Options: Optional carrying bag and carrying battery kit (7200mAh). Power

supply: Power kit (12V) to emergency vehicle, optional.

Operating system: Linux.

Dimensions: 19 x 8.2 x 10.5cm.

Weight: 1.5kg.

Warranty: 3 years.

EMERGENCY/TRANSPORT MONITOR C30 STANDARD CONFIGURATION COMEN

COMEN

Ref. BCMONI07

Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	3661040219479	HS Code	9018193010
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	26	Bare Unit / Width (cm)	8.2
Packed unit / Height (cm)	21	Bare Unit / Height (cm)	10.5
Packed unit / Length (cm)	28	Bare Unit / Length (cm)	19
Packed unit / Weight	3.2	Bare Unit / Weight (kg)	1.5

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	60
Minimum humidity	NA	Maximum humidity	93

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---