

LUXASCOPE DERMATOSCOPE LED 3.7 V (RECHARGEABLE), INCL. USB-RECHARGER


LUXAMED®

Ref. C1.526.114



Rechargeable USB

+ Graduation



18/08/2026 13:58:33 – For information purposes only, can be subject to change

LUXASCOPE DERMATOSCOPE LED 3.7 V (RECHARGEABLE), INCL. USB-RECHARGER

LUXAMED®Ref. C1.526.114

Innovative dermatoscope in pocket size, rechargeable USB.

WITH scale.

Description

- Illumination: LED-Ring with 6 LEDs in concentric structure for a highly detailed image, homogenous light
- Maintenance-free, because no lamp replacement necessary
- Removable contact plate for easy cleaning
- Stepless zoom area with 10 x magnification for a distortion-free image
- Optimized design of head
- Can be dimmed
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours by constant current drive
- Intensity of illumination: 10.000 Lux
- Colour temperature: approx. 4.200 K
- Two-component handle (aluminium / fibre-glass reinforced), also combinable with LuxaScope LED Otoscope- and Ophthalmoscope head
- Electronic ON/OFF push-button-switch
- Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive
- Incl. recharger with USB plug (adapter EU/UK/US)

**LUXASCOPE DERMATOSCOPE LED 3.7 V
(RECHARGEABLE), INCL. USB-RECHARGER****LUXAMED®****Ref. C1.526.114**

Available in

Black: ref C1.526.114

White : ref C1.526.914

Features

CE Marking	NONE	Medical device	No
Bar code	04260381963448	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.155		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Black
--------	-------