

**LUXASCOPE DERMATOSCOPE CCT LED 2.5
V - CONTACT PLATE WITH SCALE -
LUXAMED®**

Ref. C1.716.114


LUXAMED®



Technologie CCT
2 couleurs de lumière



+ Graduation
| | | | | | | | | | | | | | | | | |



18/08/2026 13:06:34 – For information purposes only, can be subject to change

**LUXASCOPE DERMATOSCOPE CCT LED 2.5
V - CONTACT PLATE WITH SCALE -
LUXAMED®**

LUXAMED®

Ref. C1.716.114

LUXASCOPE DERMATOSCOPE CCT LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED®

LUXAMED®**Ref. C1.716.114**

Innovative pocket-sized dermatoscope with CCT LED technology: choice of warm white or cool white light.

WITH graduated contact plate

Description

- Lighting: LED ring with 2X5 LEDs (concentric), no fibre optics, field of vision is not limited.
- Long press on the switch to change from warm white to cool white
- Removable contact plate for easy cleaning
- Maintenance-free, as no bulb change is required.
- Continuous zoom with X10 magnification for a distortion-free image
- Optimised head design
- Automatic switch-off after 3 minutes
- LED life: approx. 100,000 hours at constant current
- Illuminance: approx. 10,000 Lux
- Colour temperature: ring 1: approx. 3000 K, ring 2: approx. 6000 K
- Handle made of two materials: aluminium and reinforced fibreglass, compatible with the head of the Luxascope Auris CCT LED otoscope
- ON/OFF push-button switch
- Power supply: 2x AAA alkaline batteries, included

Available in :

- Black : C1.716.114
- White: C1.716.914

LUXASCOPE DERMATOSCOPE CCT LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED®

LUXAMED®
Ref. C1.716.114

Features

CE Marking	NONE	Medical device	No
Bar code	04260381963622	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.17		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		