


Ref. GLU030



CE
0537



18/08/2026 13:27:56 – For information purposes only, can be subject to change

=====

Ref. GLU030



Description

=====

Features

CE Marking	CE0537	Medical device	DMDIV (diagnostic In Vitro)
Bar code	4719869705938	HS Code	90183190

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	4	Packed unit / Height (cm)	4
Packed unit / Length (cm)	7	Packed unit / Weight	0.035

SPECIFIC STORAGE CONDITIONS

Minimum temperature	4	Max temperature	30
---------------------	---	-----------------	----

OTHER FEATURES