

THERMOMETER TEMPO 10 FLEX GREEN

Ref. SP1643/



CE
0123



THERMOMETER TEMPO 10 FLEX GREEN

Ref. SP1643/



Flexible digital thermometer, designed for accurate and comfortable measurement.

Description

The TEMPO 10 FLEX thermometer is a thermometer indicating body temperature in 10 seconds and has a flexible mouthpiece•

- Fever alarm•
- Automatic shutdown•
- Water resistant and flexible mouthpiece•
- Last measurement memory•
- Measurement in 10 seconds•
- 2 colors available green and gray
- Comes with a storage case and button battery (1.5V)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446022877	HS Code	9018191000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	6	Packed unit / Height (cm)	2
Packed unit / Length (cm)	15	Packed unit / Weight	0.052
Bare Unit / Weight (kg)	0.01		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---