

DISPLAY OF 12 THERMOMETRES TEMPO 10 FLEX GREY

Ref. SP1648P



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

DISPLAY OF 12 THERMOMETRES TEMPO 10 FLEX GREY

Ref. SP1648P



Display of 12 Thermometres Tempo 10 Flex grey

Description

Display of 12 thermometers TEMPO 10 Flex gray The thermometer TEMPO 10

FLEX is a thermometer indicating the body temperature in 10 seconds and having a flexible

- mouthpiece• Fever alarm•
- Automatic shutdown•
- Water-resistant and flexible tip•
- Last measurement memory
- Measurement in 10 seconds•
- 2 colors available green and gray• Comes with a storage case and button battery (1.5V)
- Comes with: Storage case, button battery (1.5V).

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041625	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	16
Packed unit / Length (cm)	11	Packed unit / Weight	0.66
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		