

DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO EASY GREY

Ref. SP3066P



CE
0123



DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO EASY GREY



Ref. SP3066P

The Tempo Easy allows non-contact temperature measurement through infrared rays.

Description

Display of 6 gray Tempo Easy thermometers.

The TEMPO EASY thermometer is an accurate thermometer with 3 functions for professional or home use. Non-contact technology (distance 3-5 cm) through infrared rays, thus avoiding cross-contamination. •

- 3 measurement modes: body, surface and ambient•
- Instant measurement•
- Quickly obtained temperature with audible alarm•
- 32 memories• 2
- colors available, gray or blue
- Comes with 2 LR06 batteries (AA)

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446012038	HS Code	90251920
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	20	Packed unit / Height (cm)	18
Packed unit / Length (cm)	22	Packed unit / Weight	1.92
Bare Unit / Weight (kg)	0.45		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		