

TEMP DUO II SPENGLER FOREHEAD AND AURICULAR THERMOMETER

Ref. SP3105



CE
0123



18/08/2026 12:42:20 – For information purposes only, can be subject to change

TEMP DUO II SPENGLER FOREHEAD AND AURICULAR THERMOMETER

Ref. SP3105



Temporal or ear, the measurements are instantaneous and visible on a backlit screen.

Description

- Measuring range: 32°C - 43°C

Equipped with a backlit screen, TEMPO DUO II turns red when the temperature is above normal. SPENGLER performance at home.

2 buttons allow you to take a quick and reliable measurement, reducing mishandling to zero. Change modes without any worries, calibration is done automatically.

With its small size and low weight, TEMPO DUO II slips into a bag for a trip. His 20 memoirs make it possible to follow his results over time. • Screen resolution: 0.1°C • Unit of measurement: C° or

F° • Memory: 20 measurements • Dimensions: 150 x

50 x 38.5 mm •

Weight: 75g with batteries •

Normal conditions of use: •

Ambient temperature: 16°C – 35°C

- Relative humidity: 85% or less

TEMP DUO II SPENGLER FOREHEAD AND AURICULAR THERMOMETER



Ref. SP3105

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041786	HS Code	90251920
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	8	Packed unit / Height (cm)	5
Packed unit / Length (cm)	18	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.06		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--