

NON-CONTACT THERMOMETER TEMPO PRO SPENGLER

Ref. SP3110



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



18/08/2026 12:16:56 – For information purposes only, can be subject to change

NON-CONTACT THERMOMETER TEMPO PRO SPENGLER

Ref. SP3110



Infrared (IR), non-contact, clinical thermometer, digital, handheld, range 34 # 43 °C

Description

At a distance of a few centimetres from the temple, its display instantly shows the temperature.

Without contact to avoid cross-contamination, Tempo Pro is easy to disinfect and always stays clean.

TECHNICAL CHARACTERISTICS

- Measurement : Infrared
- Distance de mesure : - 5 cm
- Unit of measurement : Celsius or Fahrenheit
- Measuring range : 32,0°C – 43,0°C / 89,6°F - 109,4°F
- Resolution : 0,1°C / F°
- Forehead accuracy : 0,2°C / 0,4F°
- Memories : 30 memories
- Measurement time: Less than 1 second
- Screen : Backlit LCD screen
- Auto power off
- 2 x 1,5V AAA batteries (included)

NON-CONTACT THERMOMETER TEMPO PRO SPENGLER

Ref. SP3110



Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041823	HS Code	9025190090
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	16	Bare Unit / Length (cm)	13
Packed unit / Weight	0.2	Bare Unit / Weight (kg)	0.13
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		