

TEMPO MAX HYPOTHERMIC DIGITAL THERMOMETER SPENGLER

Ref. THE002



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

Description

18/08/2026 13:29:31 – For information purposes only, can be subject to change

TEMPO MAX HYPOTHERMIC DIGITAL THERMOMETER SPENGLER

Ref. THE002

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6970392219271	HS Code	90251920
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	14.3
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.02

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES