

WALL UNIT MODULE WITH 1 HANDLE ONLY WITHOUT INSTRUMENT HEAD

LUXAMED®

Ref. W9.316.114



18/08/2026 13:27:30 – For information purposes only, can be subject to change

WALL UNIT MODULE WITH 1 HANDLE ONLY WITHOUT INSTRUMENT HEAD

LUXAMED®

Ref. W9.316.114

Wall Unit module with 1 handle only without instrument head

Description

The wall station is equipped to provide the required diagnostic aid, which is useful in daily life in a hospital, clinic, or private practice.

- Power supply on mains 220V / 110V : 1 power cable + 1 EU/GB/US adapter included
- Automatic shut-off after 3 min of use
- Possibility to vary the brightness with a simple long press

Available in :

- Wall station with otoscope, ophthalmoscope and speculum dispenser*
- Wall station with otoscope and speculum dispenser*

*Speculums not included with the dispenser

The device is sold with the handle only, and must be completed with an otoscope, dermatoscope or ophthalmoscope head.

WALL UNIT MODULE WITH 1 HANDLE ONLY WITHOUT INSTRUMENT HEAD

LUXAMED®Ref. W9.316.114

Features

CE Marking	NONE	Medical device	No
Bar code	4260381965022	HS Code	940549400
Country	Germany		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES