

DUAL PULSE STETHOSCOPE / BLACK

Ref. 507421



21/08/2026 00:10:34 – For information purposes only, can be subject to change

DUAL PULSE STETHOSCOPE / BLACK

Ref. 507421



Description

Pearly tube model. For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaural with built-in spring.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009205	HS Code	90181990
Warranty (months)	36		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.237

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Azure	Packed unit / Quantity	1
--------	-------	------------------------	---