

2.5V LARYNGOSCOPE HANDLE SMALL LUXAMED®

Ref. E1.416.012



2.5V LARYNGOSCOPE HANDLE SMALL LUXAMED®

Ref. E1.416.012

LUXAMED® laryngoscopes: Quality for Every Application!

Description

Laryngoscopes are used to perform intubations during surgical procedures, respiratory emergencies, or resuscitation. They are also used for therapeutic purposes to perform laryngoscopies.

The handle of the LUXAMED® laryngoscope is its key component. The quality of its lighting and its durability make it an outstanding instrument.

They are easy to clean because they are smooth and feature horizontal grooves for a better grip.

They are available with Macintosh, Miller, and McCoy blades in various sizes.

F.O. Laryngoscope, Small Handle, 2.5V LED

- Handle size: small
- Diameter: 20 mm
- Length: 142 mm
- 2.5V
- LED technology
- Stainless steel handle
- Flat surface and horizontal grooves for improved hygiene and grip.
- Quick and easy bulb replacement
- Autoclavable up to 134° (without bulb or batteries)
- Handle complies with ISO 7376 (eco-friendly standard)
- 2 alkaline batteries (AA/LR6) included

2.5V LARYNGOSCOPE HANDLE SMALL LUXAMED®

Ref. E1.416.012

Features

CE Marking	CE	Medical device	Class I
Bar code	4260381960331	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	45
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES