

## IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073\_RC



## IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073\_RC



=====

### Description

- Stethoscope with double metal horn adult (Ø: 45/30 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in single pavilion for adults

Available in Blue (DH7073\_B), Black (DH7073\_N), Red (DH7073\_R) or Pink (DH7073\_RC) for the paediatric version (Ø: 35/30 mm)

### Features

| CE Marking                         | CE            | Medical device            | Class I    |
|------------------------------------|---------------|---------------------------|------------|
| Bar code                           | 3661041000106 | HS Code                   | 9018908400 |
| Warranty (months)                  | 12            | Country                   | China      |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |            |
| Packed unit / Width (cm)           | 12            | Packed unit / Height (cm) | 4          |
| Packed unit / Length (cm)          | 25            | Packed unit / Weight      | 0.15       |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |            |
| Minimum temperature                | -20           | Max temperature           | 70         |
| Minimum humidity                   | 20            | Maximum humidity          | 85         |
| <b>OTHER FEATURES</b>              |               |                           |            |
| Colors                             | Pink          | Packed unit / Quantity    | 1          |