

## PULSE STETHOSCOPE / RED

Ref. 507522



CE



03/09/2026 23:46:47 – For information purposes only, can be subject to change

## PULSE STETHOSCOPE / RED

Ref. 507522



### Description

For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaura with built-in spring.

### Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3700446009281 | HS Code        | 90181990 |
| Warranty (months) | 36            |                |          |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |      |
|---------------------------|----|---------------------------|------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc   |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.22 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|        |           |                        |   |
|--------|-----------|------------------------|---|
| Colors | Navy blue | Packed unit / Quantity | 1 |
|--------|-----------|------------------------|---|