

STETHOSCOPE DUAL PULSE II, PINK

Ref. 507710



CE



Description

03/09/2026 23:30:08 – For information purposes only, can be subject to change

STETHOSCOPE DUAL PULSE II, PINK

Ref. 507710



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029913	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.226
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Blueberry	Packed unit / Quantity	1