

LUXASCOPE DERMATOSCOPE LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED

®

LUXAMED®

Ref. C1.516.944



+ Graduation



06/09/2026 21:08:36 – For information purposes only, can be subject to change

LUXASCOPE DERMATOSCOPE LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED ®

LUXAMED®

Ref. C1.516.944

Innovative pocket-sized dermatoscope with graduated contact plate.

Available in a range of colours.

Description

- Lighting: LED ring with 6 (concentric) LEDs, for a highly detailed image and even lighting.
- Maintenance-free, as no bulb change is required.
- Removable contact plate for easy cleaning
- Continuous zoom with X10 magnification for a distortion-free image
- Optimised head design
- Brightness can be dimmed to avoid glare
- Automatic switch-off after 3 minutes
- LED life: approx. 100,000 hours at constant current
- Illuminance: approx. 10,000 Lux
- Colour temperature: approx. 4,200 K
- Handle made of two materials: aluminium and reinforced fibreglass, compatible with the head of the Luxascope CCT otoscope and ophthalmoscope
- ON/OFF push-button switch
- Power supply: 2x AAA alkaline batteries, included

Available in :

- White: C1.516.914
- Light blue: C1.516.934

 06/09/2026 21:08:36 – For information purposes only, can be subject to change

LUXASCOPE DERMATOSCOPE LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED ®

LUXAMED®

Ref. C1.516.944

- Yellow: C1.516.944
- Fuchsia: C1.516.974
- Yellow : C1.516.944
- Kids : C1.916.914
- Black : C1.516.114
- Green : C1.516.954
- Purple : C1.516.964

Features

CE Marking	NONE	Medical device	No
Bar code	04260381963745	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	YELLOW
--------	--------