

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

Ref. M1.121.914

LUXAMED®



CE



06/09/2026 20:24:25 – For information purposes only, can be subject to change

REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

LUXAMED®Ref. M1.121.914

Reflex Hammer according to Buck | small

Description

- Reflex hammer made of aluminium
- Optimal weight distribution
- Powder coating
- With neurological sensor tip
- Length: 18 cm | weight: 90 g
- Without brush due to hygienic reasons

Available in :

White: M1.121.914

Blue :M1.121.812

Pink :M1.121.612

Yellow: M1.121.513

Kid : M1.131.914

Black M1.121.112

Green: M1.121.412

Purple: M1.121.312

**REFLEX HAMMER ACCORDING TO BUCK -
LUXAMED****LUXAMED®**Ref. M1.121.914**Features**

CE Marking	CE	Medical device	Class I
Bar code	04260381960270	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	White
--------	-------