

ULTRASOUND GEL TRANSPARENT - BOTTLE OF 1 L

Ref. CGGEL03



CE



ULTRASOUND GEL TRANSPARENT - BOTTLE OF 1 L

Ref. [CGGEL03](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	8023279332759	HS Code	30067000
Country	Turkey		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	26	Bare Unit / Height (cm)	6.5
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	27
Packed unit / Weight	1.016	Bare Unit / Weight (kg)	1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		