


SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE ADULT FINGER CLIP FOR MP3



Ref. 222260



CE
0537



10/09/2026 21:58:24 – For information purposes only, can be subject to change

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE ADULT FINGER CLIP FOR MP3



Ref. 222260

Description

SpO2 sensor compatible NELCOR OXIMAX, rigid adult for oximeter MASTER PALM 3

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030131	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.105
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

10/09/2026 21:58:24 – For information purposes only, can be subject to change