

## REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

Ref. M1.121.513

**LUXAMED®**



CE



## REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

**LUXAMED®**Ref. M1.121.513

Reflex Hammer according to Buck | small

### Description

- Reflex hammer made of aluminium
- Optimal weight distribution
- Powder coating
- With neurological sensor tip
- Length: 18 cm | weight: 90 g
- Without brush due to hygienic reasons

Available in :

White: M1.121.914

Blue :M1.121.812

Pink :M1.121.612

Yellow: M1.121.513

Kid : M1.131.914

Black M1.121.112

Green: M1.121.412

Purple: M1.121.312

# REFLEX HAMMER ACCORDING TO BUCK - LUXAMED

**LUXAMED®**
Ref. M1.121.513

## Features

|                                    |                |                           |          |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking                         | CE             | Medical device            | Class I  |
| Bar code                           | 04260381960263 | HS Code                   | 90189084 |
| Warranty (months)                  | 24             | Country                   | Germany  |
| <b>DIMENSIONS AND WEIGHT</b>       |                |                           |          |
| Packed unit / Width (cm)           | nc             | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc             | Packed unit / Weight      | nc       |
| <b>SPECIFIC STORAGE CONDITIONS</b> |                |                           |          |
| Minimum temperature                | -10            | Max temperature           | 40       |
| Minimum humidity                   | 10             | Maximum humidity          | 90       |
| <b>OTHER FEATURES</b>              |                |                           |          |
| Colors                             | YELLOW         |                           |          |