

# THE SINGLE USE PRESSURE CUFF EASY FUSE DDM

**Ref. M13000**



10/09/2026 21:57:51 – For information purposes only, can be subject to change

# THE SINGLE USE PRESSURE CUFF EASY FUSE DDM

Ref. M13000

The single use pressure cuff made in France !

## Description

The pressure infusion cuffs are medical devices designed for daily use in operating theaters, Intensive Care units or emergency transport vehicles. They can be used to manually accelerate infusion of blood and solutions when gravity alone cannot achieve the desired flow rates.

Easyfuse pressure cuffs are single use and compatible with different sizes of blood/solution bags: 500cc, 1000cc and 3000cc.

### Optimized monitoring

- **Single use**, no cleaning and no disinfection
- blood/solution bag clearly visible
- pressure indicator with automatic overpressure regulator
- Controlled decompression facilitated by 3 ways valve
- Résistant as made of TPU

All our pressure infusion cuffs are made in France in our Castelculier (Agen) manufacturing site..

## THE SINGLE USE PRESSURE CUFF EASY FUSE DDM

Ref. M13000

### Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3700468203872 | HS Code        | 98109084 |
| Warranty (months) | 12            | Country        | France   |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |      |
|---------------------------|----|---------------------------|------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc   |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.21 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES