

BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S

Ref. M05-320-250



BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S



Ref. M05-320-250

=====

Description

- Handle with tip on the screw
- Unscrewable head 2 positions
- Length 25 cm
- Available in blue (M05-320-260), black (M05-320-250) and red (M05-320-280).

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446010652	HS Code	9018908400
Warranty (months)	60	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	25
Packed unit / Weight	0.089	Bare Unit / Weight (kg)	0.089

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Black	Packed unit / Quantity	1
--------	-------	------------------------	---