

MULTY'S ELITE BAGS, BLACK

Ref. MEMULTI01_G

ELITE BAGS



13/09/2026 02:23:57 – For information purposes only, can be subject to change

MULTY'S ELITE BAGS, BLACK

Ref. MEMULTI01_G



Multi-purpose case

Description

It is designed to offer the professional the necessary support in first aid care in various professions and / or situations, such as: Physiotherapy, First Aid, Sports Activities ... - 1 padded main space, large capacity

with extractable separator- 2

padded pockets of medium capacity

- 1 flat pocket

- 2 net pockets

- 12 medium format elastics

- 8 large format elastics- Handle and shoulder strap- Color: black, blue (yellow shoulder strap)

or gray (pink shoulder strap)

- Weight: 580 gr

- Outdoor measures: 31.5 x 20 x 20 cm / 14 L

- Recommended maximum capacity: 3.5 kg

- Sold empty

- Available in black (MEMULTI01_N) and grey (MEMULTI01_G)

MULTY'S ELITE BAGS, BLACK

Ref. MEMULTI01_G



Features

CE Marking	NONE	Medical device	No
Bar code	8436564542489	HS Code	42029900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	16	Bare Unit / Width (cm)	20
Packed unit / Height (cm)	25	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	37	Bare Unit / Length (cm)	31.5
Packed unit / Weight	0.55	Bare Unit / Weight (kg)	0.58

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Gray	Packed unit / Quantity	1
--------	------	------------------------	---