

6 X BLACK BULBS WITH VALVE FOR VL CLASSIC, VL CLINIC, VL NANO, MAXI+3

Ref. 311014



Description

Classic® Black bulbs/Classic® valves

18/09/2026 01:11:27 – For information purposes only, can be subject to change

6 X BLACK BULBS WITH VALVE FOR VL CLASSIC, VL CLINIC, VL NANO, MAXI+3

Ref. 311014



Features

CE Marking	NONE	Medical device	No
Bar code	3700446026950	HS Code	9018901000
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		