

## RESCUE BACKPACK EMERGENCY, BLUE EM07.004

Ref. MEEMER05\_B

**EMERGENCY'S™**



18/09/2026 01:10:51 – For information purposes only, can be subject to change

## RESCUE BACKPACK EMERGENCY, BLUE EM07.004

**EMERGENCY'S™**

Ref. MEEMER05\_B

**Its foldable flap opening provides access to the tools necessary for emergency actions**

### Description

- Main ergonomic handle with Velcro- Large pockets on both sides
- Padded shoulder straps and backs, adjustable backpack
- Large front pocket, ready to accommodate a defibrillator
- Spacious interior for convenient and easy access to the material
- Internal compartments with mesh and transparent windows for the organization of the material
- External dimensions: 37 x 45 x 21 cm / 30L
- Material: 600D Polyester and PVC
- Weight: 1.22 kg
- Available in waterproof version
- Sold empty

Available in :

Blue : MEEMER05\_B

Red : MEEMER05\_R

# **RESCUE BACKPACK EMERGENCY, BLUE** **EM07.004**

**EMERGENCY'S™**
Ref. MEEMER05\_B

## **Features**

|                                    |               |                         |           |
|------------------------------------|---------------|-------------------------|-----------|
| CE Marking                         | NONE          | Medical device          | No        |
| Bar code                           | 8436564542663 | HS Code                 | 420299000 |
| Warranty (months)                  | 24            | Country                 | China     |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                         |           |
| Packed unit / Width (cm)           | 31            | Bare Unit / Width (cm)  | 37        |
| Packed unit / Height (cm)          | 8             | Bare Unit / Height (cm) | 21        |
| Packed unit / Length (cm)          | 48            | Bare Unit / Length (cm) | 54        |
| Packed unit / Weight               | 1.3           | Bare Unit / Weight (kg) | 1.22      |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                         |           |
| Minimum temperature                | NA            | Max temperature         | NA        |
| Minimum humidity                   | NA            | Maximum humidity        | NA        |
| <b>OTHER FEATURES</b>              |               |                         |           |
| Colors                             | Blue          | Packed unit / Quantity  | 1         |