

Ref. 332012



CE



 23/09/2026 00:57:52 – For information purposes only, can be subject to change

=====

Ref. 332012



Description

Cotton cuff with highly resistant rubber bladder.
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
 Washable at 30-40°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028480	HS Code	9018901000
Warranty (months)	6	Country	Indonesia

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.21
Bare Unit / Weight (kg)	0.21		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--