

MONOBLOCK STRAIGHT ARM CUFF WITH 1 BLADDER AND FASTENING STRAP FOR DDM PNEUMATIC TOURNIQUETS

Ref. GBDMS101

DDM 



CE



📅 29/09/2026 02:24:01 – For information purposes only, can be subject to change

MONOBLOCK STRAIGHT ARM CUFF WITH 1 BLADDER AND FASTENING STRAP FOR DDM PNEUMATIC TOURNIQUETS

Ref. GBDMS101



Monoblock straight arm cuff with one bladder and fastening strap!

Description

Monobloc PU décontaminable arm cuff are available in straight design. Made of a single monobloc unit with one bladder and fastening strap to secure the attachment.

They are available in a range of different sizes easily indentifiable by color coding of the tubing. (A : bladder width, B: bladder length, C: total length)

		A	B	C
• New born arm cuff	White tube	3.5	29	38
• Child arm cuff	Pink tube	4.5	35	48
• Small adult arm cuff	Parme tube	6	46	61
• Adult arm cuff	Sky-blue tube	8	46	67
• Large adult arm cuff	Turquoise tube	10	47	70

Monobloc arm cuffs can not be sterilized. They can be dry cleaned or immersed paying attention to maintain the tube connector out of water.

MONOBLOCK STRAIGHT ARM CUFF WITH 1 BLADDER AND FASTENING STRAP FOR DDM PNEUMATIC TOURNIQUETS



Ref. [GBDMS101](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468208877	HS Code	90189084
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES